

TCM5 FLEX monitor

Specifications

Hardware

Display

Screen: 9" color, capacitive, multi-touch TFT, 800 x 480 pixels, 16:9 ratio
 Display options: trend view, detail view
 Screen update rates: Numeric data values (SpO₂, pulse rate, tcpCO₂, tcpO₂, heating power): 1/sec
 Trend data range (depends on time axis scale): 1 pixel / 2 sec to 1 pixel / 96 sec
 Viewing angle adjustment: -10° / 0° / +8°

Operating conditions

Ambient temperature: 5-40 °C (41 - 107 °F)
 Relative humidity: up to 90 %, non-condensing
 Atmospheric pressure: 525 - 800 mmHg (700 - 1060 hPa)
 Built-in barometer:
 Range: 375-825 mmHg or 50-110 kPa
 Accuracy: ±5 mmHg or 0.67 kPa
 Power: 100-240 V, 50-60 Hz, 55-85 VA (max.)
 Rechargeable Lithium Ion battery, Type RRC2040
 Typical operating time: 4 hours per charge at 25 °C

Dimensions

TCM5 monitor

| | |
|--------|----------------------------------|
| Width | 270 mm |
| Depth | 152 mm |
| Height | 188 mm |
| Weight | 2.5 kg, 2.3 kg (without battery) |

Software

Measuring range

Transcutaneous carbon dioxide tension/tcpCO₂:
 5-200 mmHg or 0.7-26.7 kPa
 Transcutaneous oxygen tension/tcpO₂:
 0-800 mmHg or 0.0-99.9 kPa
 Oxygen saturation/SpO₂: 0-100 % (70-100 % with ±3 digits)
 Pulse rate: 25-240 bpm

Calibration

Automatic calibration
 Calibration gas: CO₂: 7.5%, O₂: 12%, N: 80.5%
 Integrated calibration chamber
 Maximum interval between two calibrations is set to 12 hours.

Patient data storage

Up to 1 year of measuring data in 1-sec data intervals
 Download of stored patient data to USB storage

Connectivity

Network connection: Ethernet 10/100 Base-T full duplex
 USB: 3x USB 2.0 (1x side, 2x rear), Type-A (USB 1.1 compatible)
 Isolated serial output: USB 2.0, Type-B (USB 1.1 compatible) & RS 232
 Isolated Nurse call output: relay contact normally open (calling, closed, not calling, opened), 60V, 2.5A
 Isolated analog output: 0-1000 mV
 Data protocol: MonLink (2.0, TCM4 series, TOSCA500, MicroGas 7650)

Interface with third-party products

Polysomnographs: Alice 5/6, Embla, Embletta Gold and other.
 Patient monitoring systems: Philips, General Electric (GE), Mindray
 Mounting system: compatible to GCX mounting system
 USB storage: Data export

Site timer

Indication of remaining measuring time
 Measuring time elapsed: clock triggers an alarm (and sensor temperature is off after 15 min)

Alarm level

OFF, 1-10 (1 is minimum)
 OFF: only available in Sleep mode

Languages

English, French, German, Czech, Danish, Dutch, Estonian, Finnish, Hungarian, Italian, Japanese, Lithuanian, Chinese, Norwegian, Polish, Portuguese, Romanian, Russian, Slovak, Spanish, Swedish, Turkish

tc Sensor 92, tc Sensor 84 and tc Sensor 54

Sensor performance (in vitro)

tc Sensor 92

Measurement principle:

Stow-Severinghaus-type tcpCO₂ sensor combined with Masimo SET® SpO₂ pulse oximetry

Conditions: sensor temperature of 43.0 °C, calibration interval: 12 h

tcpCO₂:

Accuracy: (5-200 mmHg): ± 5 mmHg

Response time: (0-90%): ≤ 70 sec

Drift: ≤ 0.5%/h

Linearity: at 1 and 10% CO₂: better than 1 mmHg or 0.13 kPa
at 33% CO₂: better than 5 mmHg or 0.67 kPa

SpO₂:

Accuracy: (70-100%): ± 3 digits

Pulse rate:

Accuracy: (25-240 bpm) ± 3 bpm

SpO₂ and pulse rate: signal averaging 2, 4, 8, 10, 12, 14 and 16 sec

Data update rate: 1 per second

tc Sensor 54

Measurement principle:

Stow-Severinghaus-type tcpCO₂ sensor

Conditions: sensor temperature of 43.0 °C, calibration interval: 12 h

tcpCO₂:

Accuracy: (5-200 mmHg): ± 5 mmHg

Response time: (0-90%): ≤ 70 sec

Drift: ≤ 0,5%/h

Linearity: at 1 and 10% CO₂: better than 1 mmHg or 0.13 kPa
at 33% CO₂: better than 3 mmHg or 0.40 kPa

tc Sensor 84

Measurement principle:

Stow-Severinghaus-type tcpCO₂ sensor combined with Clark-type tcpO₂ sensor

Conditions: sensor temperature of 43.0 °C, calibration interval: 12 h

tcpCO₂:

Accuracy: (5-200 mmHg): ± 5mmHg

Response time: (0-90%): ≤ 70 sec

Drift: ≤ 1%/h

Linearity: at 1 and 10% CO₂: better than 1 mmHg or 0.13 kPa,
at 33% CO₂: better than 5 mmHg or 0.67 kPa

tcpO₂:

Accuracy: (0-20.9 % O₂): ± 5 mmHg
(20.9-100 % O₂): ± 10 %

Response time: (0-90%): ≤ 25 sec

Drift: ≤ 1%/h

Linearity: at 0% O₂: better than 1 mmHg or 0.13 kPa
at 20.9% O₂: better than 3 mmHg or 0.4 kPa
at 50% O₂: better than 5 mmHg or 0.67 kPa
at 90% O₂: better than 25 mmHg or 3.33 kPa

Interference by anesthetic gases (in vitro)

tcpCO₂:

75% N₂O, 2% Halothane, 2% Enflurane, 2% Isoflurane,
2% Desflurane, 2% Sevoflurane: negligible

tcpO₂:

75% N₂O: <10 mmHg or 1.33 kPa

2% Halothane: approx. 200 mmHg or 26.67 kPa

2% Enflurane, 2% Isoflurane, 2% Desflurane, 2% Sevoflurane:
negligible

Sensor dimensions

Diameter: 15 mm or 0.6 in

Height: 8 mm or 0.3 in

Weight: 3 g or 0.1 oz

Sensor cable length: 3 m or 9.8 ft, shielded, flexible, polyurethane coated

Biocompatibility

Not made with natural rubber latex.

Compliance

The TCM5 monitor complies with the following standards: IEC 60601-1 (general safety), IEC 60601-1-2 (EMC), IEC 60601-1-6 (usability),

IEC 60601-1-8 (alarms), IEC 60601-2-23 (transcutaneous monitors), ISO 80601-2-61 (pulse oximeters), ISO 14971 (Risk Management), IEC 62366 (usability engineering) IEC 62304 (Software in medical devices), ISO 10993-1 (biocompatibility), EN1041 (information supplied by manufacturers), ISO 15223-1 (Symbols).

This product complies with the requirements of the Medical Device Directive 93/42/EEC.

Data subject to change without notice.

Radiometer, the Radiometer logo, ABL, AQT, TCM, RADIANCE, AQUIRE, PICO, CLINITUBES and QUALICHECK are trademarks of or used under license by Radiometer Medical ApS.

© Radiometer Medical ApS, 2700 Brønshøj, Denmark, 2020. All Rights Reserved. MAPSSS-000150 R2.

www.radiometer.com