



# Y-Knot® Flex All-Suture Anchor System

Part of CONMED's Shoulder Restoration System<sup>™</sup>. The Y-Knot<sup>®</sup> Flex for labral and capsular-based repairs helps improve surgical flexibility and access with small, all-suture anchors that deliver secure fixation with up to 80% less bone removal.\*





# Why You Should Use Y-Knot® Flex Anchors.

#### 80% Less Bone Removal

The Y-Knot 1.3mm anchor has an extremely small footprint – requiring 80% less bone removal than a 3.0mm press-fit anchor.\* In addition to providing increased flexibility when bone surface area is limited, this also preserves options if a revision surgery becomes necessary.

#### **Required Bone Removal\***

Commonly Used Labral Repair Anchors



# Required Bone Removal Comparison



#### **Closer Proximity, Stronger Repairs.**

As the number of fixation points increase, the strength of soft tissue repair has also been shown to increase<sup>1</sup>. As such, the small footprint of the Y-Knot 1.3mm anchor allows you to place multiple points of fixation in close proximity to each other which may help lead to a stronger repair.





#### 1.3MM Y-KNOT FLEX



1.8MM Y-KNOT FLEX

#### 55% More Fixation, 30% Less Creep

When deployed, Y-Knot anchors provide significant fixation advantages. The 1.3mm and 1.8 Y-Knot anchors each deliver at least 55% more fixation strength than conventional 3.0mm press fit anchors and 30% less creep under cyclic loading than competitive all-suture anchors.\*

#### Displacement (Creep) After 4000 Cycles\*

Most Commonly Used All-Suture Anchors



#### Mean Force to Failure (N)\*

Commonly Used Labral Repair Anchors



#### Minimally Invasive Entry

Slim, percutaneous drill guides require smaller skin incisions compared to conventional 3.0mm anchor instrumentation.\* The Y-Knot Flex instrumentation was engineered to help provide access with as little disruption as possible to adjacent healthy tissue.

#### SHOULDER RESTORATION SYSTEM<sup>™</sup>

NEW



1.8MM DOUBLE-LOADED (#2 SUTURE)

# The Future of Fixation

The Y-Knot<sup>®</sup> Flex Anchor System represents our next generation of arthroscopic soft tissue fixation. As we continue our commitment to developing new tools that enhance your ability to provide the highest level of minimally invasive care, Y-Knot Flex Anchors deliver a wide range of benefits to both surgeons and their patients.

## System Overview

Building on the design of our standard Y-Knot Anchor System, The Y-Knot<sup>®</sup> Flex System has features intended to help solve some of the biggest challenges that can occur during labral and capsular-based repairs:

Slim, Flexible Shaft Designed to navigate around curves and fit down a narrow percutaneous guide.

1.3mm Y-Knot Flex— Single-loaded

C.C.S

**1.8mm Y-Knot Flex**– Double-loaded (#2 Suture)

## NEW Innovative Instrumentation

We worked with surgeons to develop new and improved instrumentation that helps overcome surface area and entry limitations during these difficult repairs. Our curved and percutaneous instruments may help improve access and anchor placement.



# Y-Knot® All-Suture Anchors – Strong Purchase with 360° FormFit™ Fixation

**1.3mm Single-loaded and new 1.8mm Double-loaded (#2 Suture)** Our Y-Knot all-suture anchors feature 360° FormFit<sup>™</sup> fixation – contracting vertically while expanding laterally when deployed for strong purchase in bone with a small footprint. The 1.8mm Y-Knot Elex anchor is double-loaded with two #2 Hi-Fi<sup>®</sup> Sutures.



Y-KNOT INSERTED



Y-KNOT DEPLOYED



**FISHMOUTH GUIDE** 

SHOULDER RESTORATION SYSTEM<sup>™</sup>



**CROWN-TIPPED** GUIDE

# Instrumentation Sets

# Enhanced Access and Anchor Placement

The Y-Knot<sup>®</sup> Flex Anchor System features new instrumentation that can help to provide direct access to the surgical site and ensure ideal anchor placement. The curved, flexible instrumentation enables delivery around a curve and perpendicular access to the glenoid for posterior SLAP and inferior Bankart repairs, while our slim percutaneous guide eliminate the need for additional portals in some cases.

# Standard Instrumentation

Fishmouth and crown-tipped guide configurations are available to facilitate anchor placement on the glenoid rim or face.

### Robust and Precise

A solid shaft provides needed strength before tapering down to a narrow distal opening for precision anchor insertion.

## Drill Bits

Each ultra-small drill bit is made of flexible nitinol and has a plastic disc at the proximal end to ensure proper drilling depth.





## NEW Curved Instrumentation

This reusable guide has a distal curve to provide improved, more perpendicular access to the repair site, especially for posterior SLAP and inferior Bankart repairs.

## Percutaneous Delivery

The slim Y-Knot<sup>®</sup> Flex percutaneous T-guides can be introduced into the joint in a single step using the trocar stylet. They can alternatively be inserted using a traditional approach with a spinal needle, guide wire, and cannulated stylet.

#### **Direct Access**

The percutaneous T-Guides can be inserted in a single step or via a conventional percutaneous approach. SHOULDER RESTORATION SYSTEM™



PERCUTANEOUS SET



SHOULDER RESTORATION SYSTEM™

# Ordering Information

Description	Catalog Number
Anchors	
Y-Knot® Flex All-Suture Anchor, 1.3mm – One strand of #2 Hi-Fi® (1 white-blue)	Y1301
Y-Knot Flex All-Suture Anchor, 1.8mm – Two strands of #2 Hi-Fi® (1 white-blue, 1 white-black)	Y1802
Y-Knot Flex All-Suture Anchor, 1.8mm – Two strands of #2 Hi-Fi $^{\odot}$ (1 blue, 1 white-black)	Y1802A
Drill Bits	
1.3mm Disposable Drill Bit	Y13D
1.3mm Disposable Drill Bit, for Hard Bone	Y13DHB
1.8mm Disposable Drill Bit	Y18D
Standard Instrumentation	
Drill Guide, Fishmouth, 1.3mm	Y-G001
Drill Guide, Crown, 1.3mm	Y-G002
Blunt Obturator	Y-G003
Sharp Trocar	Y-G004
Drill Guide, Fishmouth, 1.8mm	Y-G005
Drill Guide, Crown, 1.8mm	Y-G006
Curved Instrumentation	
Curved Drill Guide, 1.3mm	Y-CG00
Curved Guide, 1.8mm	Y-CG01
Disposable Flexible Obturator, Blunt	Y-OBT1
Disposable Flexible Trocar, Sharp	Y-OBT2
Percutaneous Instrumentation	
1.3mm Disposable Percutaneous Pack with T-Guide	Y-PERC13
1.8mm Disposable Percutaneous Pack with T-Guide	Y-PERC18
Disposable Percutaneous Pack for Y-G00x Slim Reusable Guides	Y-PERC
Instrument Tray	C6178

Local 727-392-6464ConMed.comToll Free 1-800-237-0169customer\_service@conmed.com