



MacroLyte®

SureFit™

Dispersive Electrodes

SureFit™ Dual Dispersive Electrodes

Exceptional quality and value.

- For Patients from 4.4 lbs (2.2 kg)
- Pull tab to help nurses remove the pad from the patient
- Adhesive border helps to prevent fluid intrusion
- Available with or without cable

SureFit™



Catalog #	Product Description	Cable Length	Packaging
410-2000	Dual Dispersive Electrode (for patients > 4.4 lbs.), Single use, Pre-attached 10' (3.05m) Cable	10' Cable	100/case, 25/box, 1/pkg
410-2100	Dual Dispersive Electrode (for patients > 4.4 lbs.), Single use, Pre-attached 15' (4.575m) Cable	15' Cable	100/case, 25/box, 1/pkg
410-2200	Dual Dispersive Electrode (for patients > 4.4 lbs.), Single use, without cable, resealable 5-pack pouch (use reusable cable # T229)	No Cable	100/case, 5/pkg
410-2400	Dual Dispersive Electrode (for patients > 4.4 lbs.), Single use, without cable (use reusable cable # T229)	No Cable	100/case, 1/pkg

MacroLyte® Dispersive Electrodes

Cost-effective alternative with an assortment of sizes available.

- For a variety of patients from <5 kg (11 lbs) to >15 kg (33 lbs)
- Adhesive hydrogel
- Thin foam material adds in conformability
- Scalloped edges

MacroLyte®



Catalog #	Product Description	Cable Length	Packaging
400-2100	Adult Single Dispersive Electrode (for patients >15 kg), Single use, Pre-attached 10' (3.05m) Cable	10' Cable	100/case, 25/box, 1/pkg
400-2349	Adult Single Dispersive Electrode (for patients >15 kg), Single use, without Cable (use reusable cable # T229)	No Cable	100/case, 50/box, 1/pkg
400-2700	Adult Single Dispersive Electrode (for patients >15 kg), Single use, without cable, resealable 5-pack pouch (use reusable cable # T229)	No Cable	100/case, 20/box, 5/pkg
425-2200	Pediatric Single Dispersive Electrode (for patients 5-15 kg), Single use, Pre-attached 10' (3.05m) Cable	10' Cable	100/case, 25/box, 1/pkg
425-2508	Pediatric Dual Dispersive Electrode (for patients 5-15 kg), Single use, Pre-attached 10' (3.05m) Cable	10' Cable	100/case, 25/box, 1/pkg
425-2725	Pediatric Single Dispersive Electrode (for patients 5-15 kg), Single use, without Cable (use reusable cable # T229)	No Cable	100/case, 50/box, 1/pkg
425-2801	Pediatric Dual Dispersive Electrode (for patients 5-15 kg), Single use, without Cable (use reusable cable # T229)	No Cable	100/case, 25/box, 1/pkg
440-2400	Premie Single Dispersive Electrode (for patients < 5 kg), Single use, Pre-attached 10' (3.05m) Cable	10' Cable	40/case, 10/box, 1/pkg
450-2300	Infant Single Dispersive Electrode (for patients 5-15 kg), Single use, Pre-attached 10' (3.05m) Cable	10' Cable	60/case, 15/box, 1/pkg

ConMed Corporation offers a Hold Harmless Agreement. Please contact ConMed Customer Service for more information.

CONMED CORPORATION PRODUCT AREAS:

ARTHROSCOPY • ELECTROSURGERY • ENDOSCOPY • ENDOSURGERY • GASTROENTEROLOGY • INTEGRATED SYSTEMS • PATIENT CARE • POWERED INSTRUMENTS • PULMONOLOGY



525 French Road
 Utica, NY 13502
 (315) 797-8375
 Customer Service: 1-800-448-6506
 FAX: 1-800-438-3051
 International Sales: +1(315) 797-8375 • FAX: +1(315) 735-6235
 email: info@mail.conmed.com
 www.conmed.com