




BATRIK INSTRUMENT PROTECTORS

Our Instrument Protectors protect the tips and edges of instruments to prevent damage during handling and sterilization. Available in a wide array of color coded sizes for almost every application, including osteotome, retractor, double, and perforated.

- Steam and gas permeable
- Perforated protectors reduce vacuum seal and moisture build-up
- Validated for pre-vacuum sterilization
- Multiple color-coded sizes
- Latex Free
- Non-Sterile




PRODUCT SELECTION

	ITEM # REGULAR	ITEM # PERFORATED	COLOR	DIMENSIONS
 RC-07  PC-07	REGULAR & PERFORATED / CLEAR			
	RC-01	PC-01	Clear	0.062" x 0.75" (1.57 mm x 19 mm)
	RC-02	PC-02	Clear	0.075" x 0.75" (2 mm x 19 mm)
	RC-03	PC-03	Clear	0.109" x 0.75" (2.77 mm x 19 mm)
	RC-04	PC-04	Clear	0.125" x 1" (3.18 mm x 25.4 mm)
	RC-05	PC-05	Clear	0.187" x 1" (4.75 mm x 25.4 mm)
	RC-06	PC-06	Clear	0.406" x 0.75" (10.3 mm x 19 mm)
	RC-07	PC-07	Clear	0.062" x 0.375" x 1" (1.57 mm x 9.53 mm x 25.4 mm)
	RC-08	PC-08	Clear	0.078" x 0.625" x 1" (2 mm x 15.875 mm x 25.4 mm)
	RC-09	PC-09	Clear	0.125" x 1" x 1" (3.18 mm x 25.4 mm x 25.4 mm)
	RC-10	PC-10	Clear	Assorted
PACKAGING: Bag of 100				
 RT-07  PT-07	REGULAR & PERFORATED / TINTED			
	RT-01	PT-01	Tinted White	0.062" x 0.75" (1.57 mm x 19 mm)
	RT-02	PT-02	Tinted Blue	0.075" x 0.75" (2 mm x 19 mm)
	RT-03	PT-03	Tinted Green	0.109" x 0.75" (2.77 mm x 19 mm)
	RT-04	PT-04	Tinted Red	0.125" x 1" (3.18 mm x 25.4 mm)
	RT-05	PT-05	Tinted Yellow	0.187" x 1" (4.75 mm x 25.4 mm)
	RT-06	PT-06	Tinted Maroon	0.406" x 0.75" (10.3 mm x 19 mm)
	RT-07	PT-07	Tinted Orange	0.062" x 0.375" x 1" (1.57 mm x 9.53 mm x 25.4 mm)
	RT-08	PT-08	Tinted Brown	0.078" x 0.625" x 1" (2 mm x 15.875 mm x 25.4 mm)
	RT-09	PT-09	Tinted Black	0.125" x 1" x 1" (3.18 mm x 25.4 mm x 25.4 mm)
	RT-10	PT-10	Tinted Assorted	Assorted
PACKAGING: Bag of 100				
 RS-07  PS-07	REGULAR & PERFORATED / SOLID			
	RS-01	PS-01	Solid White	0.062" x 0.75" (1.57 mm x 19 mm)
	RS-02	PS-02	Solid Blue	0.075" x 0.75" (2 mm x 19 mm)
	RS-03	PS-03	Solid Green	0.109" x 0.75" (2.77 mm x 19 mm)
	RS-04	PS-04	Solid Red	0.125" x 1" (3.18 mm x 25.4 mm)
	RS-05	PS-05	Solid Yellow	0.187" x 1" (4.75 mm x 25.4 mm)
	RS-06	PS-06	Solid Maroon	0.406" x 0.75" (10.3 mm x 19 mm)
	RS-07	PS-07	Solid Orange	0.062" x 0.375" x 1" (1.57 mm x 9.53 mm x 25.4 mm)
	RS-08	PS-08	Solid Brown	0.078" x 0.625" x 1" (2 mm x 15.875 mm x 25.4 mm)
	RS-09	PS-09	Solid Black	0.125" x 1" x 1" (3.18 mm x 25.4 mm x 25.4 mm)
	RS-10	PS-10	Solid Assorted	Assorted
PACKAGING: Bag of 100				

ITEM #	DIMENSIONS
OP-11	0.25" x 0.75" (6.35 mm x 19.05 mm)
OP-12	0.375" x 1" (9.53 mm x 25.4 mm)
OP-13	0.5" x 1" (12.7 mm x 25.4 mm)
OP-14	0.75" x 1" (19.05 mm x 25.4 mm)
OP-15	1" x 1" (25.4 mm x 25.4 mm)
OP-16	1.25" x 1.25" (31.75 mm x 31.75 mm)
OP-17	1.5" x 1.5" (38.1 mm x 38.1 mm)
OP-18	1.75" x 1.25" (44.45 mm x 31.75 mm)
OP-20	Assorted
PACKAGING: Bag of 50	

OSTEOTOME PROTECTORS




ITEM #	DIMENSIONS	PACKAGING
RGS-01	1.5" (38 mm)	55 / Bag
RGS-02	2" (51 mm)	45 / Bag
RGS-03	3" (76.2 mm)	40 / Bag

SINGLE-TIP RETRACTOR PROTECTORS



ITEM #	COLOR	SIZE
RGM-01	Pink	Small
RGM-02	Clear	Small
RGM-03	Green	Medium
RGM-04	Clear	Medium
PACKAGING: Bag of 50		

WEITLANER RETRACTOR PROTECTORS



ITEM #	COLOR	DIMENSIONS
DG-01	Clear	0.0625" x 0.75" (1.59 mm x 19.05 mm)
DG-02	Blue	0.078" x 0.75" (1.98 mm x 19.05 mm)
DG-03	Green	0.109" x 0.75" (2.78 mm x 19.05 mm)
DG-04	Red	0.125" x 1" (3.18 mm x 25.4 mm)
DG-05	Yellow	0.187" x 1" (4.76 mm x 25.4 mm)
DG-07	Orange	0.062" x 0.375" x 1" (1.59 mm x 9.53 mm x 25.4 mm)
DG-10	Assorted	Assorted
PACKAGING: Bag of 50		

PERFORATED DOUBLE PROTECTORS

