

Lanza™ Tapered Specialty Shaver Blade

Now available in 15° pre-bent shoulder and knee, and hip configurations.

CONMED's NEW Lanza[™] Tapered Specialty Shaver Blade has a tapered distal tip, smooth cutting profile and dual-inner cutting window that is designed to provide precise, efficient resection in tight anatomical spaces and sensitive areas.

Features and Benefits

Better Access	Tapered Distal Tip – Designed to provide easy access to tight spaces, difficult-to-reach anatomy and around sensitive areas.	
Faster Cutting	Dual Inner Cutting Window – Compared to a standard blade ¹ , the dual inner cutting window provides twice as many cuts per rotation to help maximize volumetric resection.	
Reduced Particulate	TrimRite® Steel & Shaft Lubrication – Designed to reduce metal particulate generation, the shrink tubing on the inner blade shaft provides a barrier to prevent metal-on-metal contact ² , while the exterior features hardened TrimRite® stainless steel that helps maintain a sharp cutting edge and blade integrity. ³	
Smooth Cutting	Smooth Cutting Profile – Engineered to prevent snagging on soft tissue when navigating tight spaces, while the shallow exterior window and dual inner cutting window are designed to reduce clogging.	
15° Pre-bent Configurations	13cm and 19cm[*] Lengths – Now available in shoulder and knee, and hip length configurations. Designed to provide precise, efficient resection, the two 15° pre-bent options make it easier to navigate around tight anatomical spaces and sensitive structures.	
	Lanza™ can also serve as an effective soft tissue probe.	

NEW



525 French Road, Utica, NY 13502

LANZA™ TAPERED SPECIALTY BLADES



Lanza[™] Blades are compatible with CONMED's Advantage® Turbo and Ergo[™] Shaver Handpieces and our D3000 and D4000 consoles.

Ordering Information

	Description	Catalog Number
	Lanza [™] Tapered Specialty Shaver Blade	9391A
NEW	Pre-bent Lanza™, 13cm	9391M
NEW	Pre-bent Lanza™, 19cm*	HPS-CB09
	Advantage® Turbo Shaver Handpiece	
	Two-Button Advantage® Turbo Shaver	D9924
	No-Button Advantage® Turbo Shaver (Footswitch Only)	D9920
	Ergo™ Shaver Handpiece	
	Ergo™ Shaver Handpiece, No-Button	D4200
	Ergo™ Shaver Handpiece, Two-Button	D4240
	Electric Consoles and Accessories	
	D3000 Arthro Power Console	D3000
	D4000 Arthro Power Console	D4000

*Not for sale at this time.

Local 727-392-6464 Toll Free 1-800-237-0169 CONMED.com customer_service@CONMED.com