

**Infiltrationspumpe für die Liposuction
Infiltration pump for Liposuction
Bomba de infiltración para la liposucción
Pompe d'infiltration pour la liposuction
Pompa d'infiltrazione per la liposuzione**

LS-i30

INFILTRATION



© Copyright 06/2015



SURGICAL

medicon®

Medicon eG · Gänsäcker 15 · D-78532 Tuttlingen
Tel.: +49 (0) 74 62 / 20 09-0 · Fax: +49 (0) 74 62 / 20 09-50
E-Mail: sales@medicon.de · Internet: www.medicon.de · Germany

Infiltrationspumpe für die Liposuction
Infiltration pump for Liposuction
Bomba de infiltración para la liposucción
Pompe d'infiltration pour la liposuction
Pompa d'infiltrazione per la liposuzione

LS-i 30
INFILTRATION



Besonders geeignet für die Injektion der
Tumesenzlösung für die Liposuktion.

Particularly suitable for the injection of tumescent
solution during liposuction procedures.

Especialmente apto para la inyección de la solución
tumescente para la liposucción.

Particulièrement approprié pour l'injection de la
solution d'infiltration pour la liposuction.

Specialmente adatta all'iniezione della soluzione
tumescente per la liposuzione.

Spannung (umschaltbar)
Voltage (switchable)
Voltaje (conmutable)
Tension (réglable)
Tensione (commutabile)

Pumpenleistung
Pumpage
Rendimiento de la bomba
Puissance de la pompe
Potenza della pompa

Gewicht netto
Net weight
Peso neto
Poids net
Peso netto

LS-i 30 INFILTRATION	100/115/230 V / 50 – 60 Hz	1 – 15 l/h	2,4 kg
----------------------	----------------------------	------------	--------

LS-i 30 Infiltration komplett:
complete:
completo:
complet:
completo:

91.70.30



91.70.32
26 x 25 x 10 cm

1

Steuergerät 100/115/230 V / 50-60 Hz
Control unit 100/115/230 V / 50-60 Hz
Unidad de control 100/115/230 V / 50-60 Hz
Unité centrale 100/115/230 V / 50-60 Hz
Unità di comando 100/115/230 V / 50-60 Hz



91.70.31

1

ON/OFF - Pedal
ON/OFF - Pedal
Pedal - ON/OFF
Pédale - ON/OFF
Pedale - ON/OFF



91.70.23

1

Flaschenhalter
Jar holder
Porta botellas
Porte-bocaux
Portavasi

Zubehör
accessories
Accesorios
Accessoires
Accessori

91.70.03

Schlauchset, einweg, steril, P/10
Set of disposable tubes, single use, sterile, P/10
Juego de tubos, desechable, estéril, P/10
Jeu de tuyaux, usage unique, stérile, P/10
Set di tubi, mono uso, sterile, P/10