

Two-Spike



Disposable Gravity Irrigation Tubing Sets

Both the Two-Spike and Four-Spike Gravity Tubing Sets are compatible with existing scope/inflow sheaths or cannulas. No special sheaths or cannulas are required. Extra large tubing diameters maintain high flow rates under the most demanding conditions.

The Four-Spike tubing set is designed to be used for high fluid demands by reducing the need for staff to constantly replenish fluid supply.

Tubing sets are packaged sterile, single use, eight units per box. Each tubing set includes the following:

- Tubing material that is latex free, made from flexible polyvinyl chloride.
- Inflow line with flexible connection (11 ft.).
- Spike line with two or four bag spikes for irrigation supply.
- Inflow clamp within sterile field to control flow rates.

Cat. No.

Gravity Irrigation Disposable Tubing Set (Two-Spike)	C7110
Gravity Irrigation Disposable Tubing Set (Four-Spike)	C7111

For more information on this and other Linvatec products, call your local representative or Linvatec Corporation.



11311 Concept Boulevard, Largo, Florida 33773-4908
 813-392-6464
 1-800-237-0169
 International Fax: 813-397-4540
 Fax: 813-399-5256

©1998 Linvatec Corporation.
 A subsidiary of ConMed® Corporation.
 All rights reserved. CBR0017

Four-Spike

