

Stealth[™] Diamond-Like Carbon Coated Shaver Blades

CONMED's top-of-the-line performance blades, Stealth[®] arthroscopic shaver blades feature a hard, gall-resistant Diamond-Like Carbon coating to reduce metal particulate generation,^{1,2} lubrication on the inner tube to reduce metal-on-metal contact, and our patented XtraSharp[®] technology for optimum cutting performance.

Features and Benefits

| Reduced Metal Particulate Generation | Diamond-Like Carbon coating on distal cutting tip has a lower coefficient of friction than standard steel and has shown reduced metal particulate generation when compared to leading competitors.^{1.2.5} Lubrication (full length shrink coverage) on the inner tube provides a barrier designed to prevent metal-on-metal contact – further helping reduce metal particulate generation within the blade shaft. |
|--|---|
| Cleaner, More Efficient Cut | Two-piece inner cutting tip design provides a 22% tighter clearance between the inner and outer cutting window allowing for a cleaner, more efficient cut. ^{3,4} |
| Optimum Cutting Performance | Patented XtraSharp[™] electrochemical grinding method provides a sharp cutting edge for ideal cutting performance. Stealth's Rigid TrimRite[®] Cutting Tip is a hardened stainless steel utilized to maintain a sharp cutting edge and optimize cutting efficiency.⁴ |



STEALTH[™] DIAMOND-LIKE CARBON COATED SHAVER BLADES

The design features of the Stealth[™] DLC Coated Shaver Blades help to reduce metal particulate generation against some leading competitors – as shown in the data comparison chart below.

Metal Particulate Comparison Chart*



*TR17-00073, TR17-00165, TR17-00073-1. Competitor A - Blade 1: Smith & Nephew, Dyonics Bonecutter. Competitor A - Blade 2: Smith & Nephew, Incisor.



Ordering Information

| Description | Size | Length | Qty. | Catalog No. |
|---|----------|--------|-------|-------------|
| Stealth [™] Diamond-Like Carbon (DLC) Coated | d Blades | | | |
| Stealth [™] Blade, Great White® | 4.2mm | 13cm | 1 box | S9299A |
| Stealth™ Blade, UltraFRR® | 4.2mm | 13cm | 1 box | SDPS-C010 |
| Stealth™ Blade, FRR | 4.2mm | 13cm | 1 box | S9247A |
| Stealth™ Blade, UltraGator® | 4.2mm | 13cm | 1 box | SDPS-C002 |
| Stealth [™] Blade, Gator® | 4.2mm | 13cm | 1 box | S9263A |

CONMED Corporation 525 French Road, Utica, NY 13502 Local 727-392-6464 Toll Free 1-800-237-0169 CONMED.com customer_service@conmed.com