

# PlumePen<sup>®</sup> Ultra Surgical Smoke Evacuation Pencil

No matter the surgical procedure, PlumePen® Ultra is up to the challenge. An addition to the PlumePen portfolio, Ultra offers a telescopic capture port with an advanced swivel design that will delight the surgeon with superior hand movement and less wrist drag.





Suffato PLTER



PLUMEPEN® ULTRA SURGICAL SMOKE EVACUATION PENCIL

### **Features & Benefits**

#### **Telescopic Capture Port with Enclosed Blade**

• Blade moves in unison with capture port affording versatility in any surgical situation. Integrated design eliminates need for blade replacement saving time and money

#### **Cam Lock System**

 Prevents capture port movement during surgery. Audible click upon activation/ deactivation of locking mechanism

#### **Responsive, Multi-Directional Swivel**

· Unique, floating ball joint design allows unrestricted hand movement

#### Integrated Wire/Tubing Management

• Ultra-light, PVC free, flexible tubing designed to reduce wrist pull and hand fatigue



```
Telescopic Capture Port with Enclosed Blade
```

To learn more about this and other innovative products, call **1-800-448-6506** or visit **CONMED.com**.

## Ordering Information

	Description	Qty	Catalog Number
NEW	PlumePen <sup>®</sup> Ultra (Non-Stick ESU Blade) with 3/8 in (9.5 mm) x 10 ft (3 m) hose/ tubing & 7/8 in (22mm) filter port adapter	20/Case	PLPUL2020
	PlumePen <sup>®</sup> Ultra (Non-Stick ESU Blade) with 3/8 in (9.5 mm) x 15 ft (4.6 m) hose/ tubing & 7/8 in (22 mm) filter port adapter	20/Case	PLPUL2520

Distributed by: CONMED Corporation 525 French Road, Utica, NY 13502

Customer Service: 1-800-448-6506 International Sales: 1-315-797-8375 Fax: 1-800-438-3051 CONMED.com info@mail.CONMED.com