

181

L,

## MAKE SMARTER, MORE INFORMED CLINICAL DECISIONS.

## Meet the Welch Allyn Connex Spot Monitor

It's easy to use and supports comprehensive bedside monitoring—from routine vital signs acquisition to early warning scoring and interval monitoring. And we can integrate with your EMR, giving your care teams immediate access to critical information—when and where they need it.



## DESIGNED WITH PATIENT OUTCOMES IN MIND.

Customized features help clinicians identify and respond to signs of patient deterioration sooner

Vital signs and other patient health data points are sent directly to the EMR to help reduce errors and omissions

3

1

2

Clinical technologies support and produce readings you can rely on

## VISUALLY ASSESSING RESPIRATION RATE? GO DIGITAL.

Digitally spot-check a patient's respiration rate without changing your workflow. Acquire accurate, connected readings at the fingertip helping your team reduce transcription errors and detect signs of deterioration.

# SECURE YOUR NETWORK AND PATIENT HEALTH DATA.

DOD RMF APPROVED'

If data security is your top concern, we have you covered. When you equip your **Connex** Spot Monitor with **Welch Allyn** EMR connectivity solutions, we align with your network security protocols to help you protect your patients and infrastructure.

- FIPS-approved<sup>2</sup> encryption of transmitted data with device and server authentication
- Digitally signed software images and configuration files protect devices from unauthorized changes
- Single sign-on compatibility, including **Imprivata**, supports quick, secure logins
- Screen lock feature protects idle devices from unauthorized access
- Patient data is automatically deleted after it's sent to the EMR and cleared from the device

### ACCESSORIES

7000-APM	Accessory Power Management Stand
7000-MS3	Classic Mobile Stand
7000-MWS	Mobile Work Surface Stand
7000-GCX	GCX Wall Channel Mount
6000-NC	VSM 6000 Nurse Call Cable
7000-916HS	2D Barcode Scanner
7000-916HSR	2D Barcode Scanner + High Frequency RFID Reader
BATT22	CSM Lithium-Ion Battery
BATT99	APM Stand Lithium-Ion Battery
4500-35	Blood Pressure Hose with <b>FlexiPort</b> (10 ft)
7000-33	Neonatal Blood Pressure Hose (10 ft)
S1-CSM-PR0-PS	CSM, SmartCare Service Protection Program, 3-years
S1-CSM-PROPLUS-PS	CSM, SmartCare Protection+ Program, 3-years

The **Connex** Spot Monitor is compatible with a wide range of cleaning products to help you reduce the risk of cross-contamination.



### ACCESSORY POWER MANAGEMENT STAND

### 7000-APM

Extend your battery life to up to 17 hours of use, helping minimize downtime due to re-charging. Plus, the integrated work surface and built-in accessory/cleaning storage were designed to enhance your productivity.

## Want to learn more? Contact your Hillrom account representative today for a demo or visit hillrom.com/csm.

### References

1. Medical Device and Equipment Risk Assessment. Data on File.

2. Validation OpenSSL FIPS Object Module v2.0. Data on File.

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

Baxter, Connex, FlexiPort, Hillrom, Welch Allyn are trademarks of Baxter International Inc. or its subsidiaries.

Medtronic and Microstream are registered trademarks of Medtronic, Inc. Masimo is a registered trademark of Masimo Corp.

Imprivata is a trademark of Imprivata, Inc.

US-FLC45-230022 V1 05/2023

## Baxter

The FDA issued Enforcement Policy for Non-Invasive Remote Monitoring Devices Used to Support Patient Monitoring During the Coronavirus Disease 2019 (COVID-19) Public Health Emergency. During this emergency and while the policy is in effect, FDA does not intend to object to limited modifications to the FDA-cleared indications without prior submission of a 510(K) where the modification does not create an undue risk. Hillrom does not yet have FDA 510(k) clearance on the combination use of the Connex Spot Monitor with Masimo RRp. Hillrom intends to adhere to FDA's recommendations to market CSM with Masimo RRp with appropriate testing and labeling while the policy is in effect.

This device is intended to provide recommendations that should be used in an adjunctive (supportive) manner and are not intended to be used as a primary means to make diagnosis, prevention, or treatment recommendations.

## Modifications to FDA Cleared Indications for Use (modifications are underlined)

The Connex Spot Monitors are intended to be used by clinicians and medically qualified personnel for monitoring of noninvasive blood pressure, pulse rate, noninvasive functional oxygen saturation of arteriolar hemoglobin (SpO<sub>2</sub>), and body temperature in normal and axillary modes of neonatal, pediatric, and adult patients. <u>Monitoring Respiration Rate from photophlethysmogram (RRp) is indicated for adult and pediatric patients greater than 2 years old.</u>

The most likely locations for patients to be monitored are general medical or surgical floors and general hospital and alternate care environments. The product is available for sale only upon the order of a physician or licensed health care professional.

### WELCH ALLYN CONNEX SPOT MONITOR SPECIAL INFORMATION FOR US ONLY

#### **Device Performance**

Validation of the integration of Masimo RRp technology into the CSM device was completed through software verification testing and design validation of the RRp parameter in the device user interface and IFU. The CSM device has been tested and shown to comply with IEC 60601-1 Edition 3.1 and IEC 60601-2-2 4th Edition. A Risk assessment has been performed according to ISO 14971. Any identified hazards have been found to be acceptable.

### **Potential Risks**

See the Instructions for Use included on the enclosed CD for a complete list of Warnings and Cautions.

For further information on the Hillrom Welch Allyn Connex Spot Monitor, including the Instructions For Use, please visit hillrom.com.