

ProTwin

Professional spray system – Trolley wash unit

APPLICATION

Designed for reprocessing Hospital S/S equipment/furniture e.g.

Stainless steel trolleys, close transport cabinets etc.

Use in Central Sterile Processing/CSSD and the OR/point-of-use Sterile Processing/TSSU located within a central surgery department.

The trolley wash unit works equally well at healthcare centers or smaller clinics, where an effective and hygienic trolley wash unit is required.

Time-saving, hygienic manual trolley washer with wall-mounted spray-gun unit, including holder for detergents.

ProTwin is a cleaning and sanitizing station that automatically mixes and dispenses concentrated chemicals. With a clean, professional appearance it comes with a built-in hose storage.



Performance

- ✦ Pressure: Min 1 bar Max6 bar (90 psi)
- ✦ Flow-rate (Rinse) 2.7 gpm (Prod) 3.5 gpm
- ✦ ProTwin foam kit
- ✦ Utilizes the ProTwin's single selector for detergent and sanitizing functions
- ✦ ProTwin's venturi maintain a proper and consistent chemical action/consumption
- ✦ Compressed medical air facilitate drying of items (to be installed)
- ✦ Multi-layer hose, improves flexibility and chemical resistant,

Adjustable spray pattern & thick foam all in one



Features

- ✦ Adjustable spray pattern
- ✦ Separate selectable chemical position eliminate risk of chemical crossover
- ✦ Adjustable accessories – easy to use
- ✦ Venturi ensure accurate dosing, easy to clean and be replaced
- ✦ Intuitive ergonomic design for user friendly operation
- ✦ 3 functions spray/foam and rinse up to 60°C (140°F) tolerance
- ✦ Multi-layer hose, improves flexibility and chemical resistant