

RZ Medizintechnik GmbH
Unter Hasslen 20 | D-78532 Tuttlingen

Phone: +49 (0)7462/9470-0

Fax: +49 (0)7462/9470-50

sales@rz-medizintechnik.com


www.rz-medizintechnik.com



RZ
Medizintechnik
GmbH

progress in surgery


Genuine TABLE MOUNTED Ring Retractor System



GENUINE RING RETRACTOR SETS


The genuine and established RZ Table Mounted Ring Retractor System provides a large variety of configurations. Every kit can be easily set up and adjusted multi-directional even for complex requirements thanks to flexible ratchet mechanisms. Any retractor can be applied to make the system an individual set.

- Large assortment of retractor blades in different sizes
- Expandable basic components to assemble a complex system
- Simple assembly and modifications before and during surgery
- Cost-effective due to reduced operation time
- All parts fully autoclavable




RETRACTOR KIT

170-800-904	Expandable general kit for wide application field	Page 6 - 7
-------------	---	------------




RETRACTOR KIT 18"

170-800-900	Expandable general kit for bariatric surgeries	Page 8 - 9
-------------	--	------------




RETRACTOR KIT 24"

170-800-902	Expandable general kit for bariatric surgeries	Page 10 - 11
-------------	--	--------------




RECTAL KIT

170-800-906	Special rectal kit based on general set	Page 13
-------------	---	---------




VAGINAL KIT STANDARD / ALTERNATIVE

170-800-909 170-800-910	Vaginal Set with special angled frames	Page 14 - 15
----------------------------	--	--------------




IMA KIT

170-800-905	Internal mammary artery access	Page 16
-------------	--------------------------------	---------




SMALL INCISION KIT

170-800-907	Hernia repair, appendectomy, cholecystectomy, vascular procedures	Page 17
-------------	---	---------



SPINAL KIT

170-800-908	Open field Retractor System	Page 18 - 19
-------------	-----------------------------	--------------



STERILISATION

	Organisation and storage of all components	Page 20 - 21
--	--	--------------




ACCESSORIES

	Table Post, Ratchets, Frames and Retractors	Page 22 - 29
--	---	--------------


GENUINE RING RETRACTOR




Height- and horizontal adjustable




Safe and solid fixture to table and/or splint



Multidirectional hinges for complex set-ups



Solid and flexible footing of frame segments



Different bar designs allow standard and complex settings




Sets can be combined with every blade type




Segments height-adjustable and expandable




Large variety of ratchet mechanisms for maximum flexibility



Stable and reliable design, fully autoclavable



Storage baskets available





















RING RETRACTOR KIT

Flexible and expandable Retractor kit for a wide field of applications. Ensures fast and simple readjustment, repositioning and replacement.

170-800-904

RETRACTOR KIT





BARIATRIC RETRACTOR SET 18"

170-800-900


BARIATRIC RETRACTOR SET

Flexible and expandable Retractor kit for a wide field of applications. Ensures fast and simple readjustment, repositioning and replacement. Especially for big or obese patients.


18"




170-800-099 Page 22 1x
Table Post for fixation at OP table




170-800-102 Page 22 1x
Horizontal Bar




170-800-110 Page 23 6x
Tilt Ratchet Mechanism




170-800-104 Page 22 1x
Post Coupler




170-800-106 Page 23 1x
Ratchet, rotatable




170-800-105 Page 23 1x
Ratchet Extender




170-800-012 Page 24 1x
Frame, oval, medium




170-800-234 Page 26 1x
Malleable Retractor Blade




170-800-235 Page 26 1x
Malleable Retractor Blade




170-800-236 Page 26 1x
Malleable Retractor Blade




170-800-214 Page 28 1x
Kelly Retractor




170-800-215 Page 28 1x
Kelly Retractor



170-800-218 Page 28 1x
Kelly Retractor Blade



170-800-201 Page 27 1x
Balfour Retractor



170-800-222 Page 28 1x
Harrington Retractor

















BARIATRIC RETRACTOR SET 24"

170-800-902

BARIATRIC RETRACTOR SET


Flexible and expandable Retractor kit for a wide field of applications. Ensures fast and simple readjustment, repositioning and replacement. Especially for big or obese patients.

24"




170-800-901 18"
170-800-903 24"

460mm, included in 170-800-901 (18")



170-800-099 Page 22 1x
Table Post for fixation at OP table


610mm, included in 170-800-903 (24")




170-800-098 Page 22 1x
Table Post for fixation at OP table




170-800-105 Page 23 1x
Ratchet Extender




170-800-106 Page 23 1x
Ratchet, rotatable



170-800-236 Page 26 1x
Malleable Retractor Blade



170-800-234 Page 26 1x
Malleable Retractor Blade



170-800-235 Page 26 1x
Malleable Retractor Blade



170-800-218 Page 28 1x
Kelly Retractor Blade




170-800-201 Page 27 1x
Balfour Retractor




170-800-222 Page 28 1x
Harrington Retractor, modified

To be used with standard Ring kit including table post, horizontal bar, couples and ratchets.




170-800-906



170-800-021 Page 24 1x
Frame, round, small



170-800-311 Page 27 2x
Rectal Blade, small



170-800-312 Page 27 2x
Rectal Blade



VAGINAL RETRACTOR SET


For optimal exposure, this vaginal retractor set includes a special 20° angled and ergonomic designed frame.

170-800-909


VAGINAL RETRACTOR SET

STANDARD SET-UP

		
170-800-100 Page 22 1x Table Post for fixation at OP table	170-800-103 Page 22 1x Horizontal Flex Bar, 360° rotatable	170-800-110 Page 23 4x Tilt Ratchet Mechanism
		
170-800-104 Page 22 1x Post Coupler	170-800-051 Page 25 1x Vaginal Frame, oval, 20° angled	170-800-261 Page 29 1x Deaver Retractor, small
		
170-800-262 Page 29 1x Deaver Retractor, large	170-800-271 Page 29 2x Vaginal Retractor	170-800-272 Page 29 1x Vaginal Retractor









ALT. SET-UP 170-800-910

		
170-800-097 Page 22 1x Table Post for fixation at OP table	170-800-112 Page 23 4x Rotate- and tiltable Ratchet Mechanism	170-800-050 Page 25 1x Vaginal Frame, oval, segmented
		
170-800-271 Page 29 2x Vaginal Retractor	170-800-237 Page 26 2x Malleable Retractor Blade, with lip	

IMA RETRACTOR KIT

Ensures optimal and flexible holding and efficient dissecting of the internal mammary artery.





170-800-905



SMALL INCISION KIT

Very flexible and stabile system for retraction, hernia repair, appendectomy, cholecystectomy and some vascular procedures.

170-800-907





SPINAL RETRACTOR SET


Creates an open field for the surgeon and allows the operating staff hands-free assistance.

170-800-908


SPINAL RETRACTOR SET




170-800-100 Page 22 1x
Table Post for fixation at OP table




170-800-103 Page 22 1x
Horizontal Flex Bar




170-800-104 Page 22 1x
Post Coupler




170-800-041 Page 24 1x
Ring Segment, round, small, 2 pieces



170-800-061 Page 25 1x
Ring Segment, small, 2 pieces




170-800-062 Page 25 1x
Ring Segment, medium, 2 pieces




170-800-063 Page 25 1x
Ring Segment, large, 2 pieces




170-800-321 Page 27 2x
Scoville Hook, small




170-800-111 Page 23 6x
Self-Retaining Tilt Ratchet




170-800-110 Page 23 4x
Tilt Ratchet Mechanism




170-800-322 Page 27 2x
Scoville Hook, large




170-800-331 Page 29 2x
Meyerding Blade, small




170-800-332 Page 29 2x
Meyerding Blade, medium




170-800-333 Page 29 2x
Meyerding Blade, large



170-800-343 Page 27 2x
Rake Retractor, 3 prongs





170-800-344 Page 27 2x
Rake Retractor, 4 prongs



170-800-345 Page 27 2x
Rake Retractor, 6 prongs

STERILISATION CASES


Easy and convenient organisation of all components for sterilisation procedures and safe storage.




Safety clasp system




Practical placing of the complete set



Storage and sterilization preparation in one basket



Easy fixation of each part



STERILISATION

170-800-912 SET


Article No.	Description	Page	Pieces
170-800-113	Standard Ratchet Mechanism	23	2
170-800-110	Tilt Ratchet Mechanism	23	6
170-800-104	Post Coupler	22	1
170-800-221	Harrington Retractor	28	1
170-800-214	Kelly Retractor, 76mm	28	2
170-800-215	Kelly Retractor, 102mm	28	1
170-800-216	Kelly Retractor, 127mm	28	1
170-800-217	Kelly Retractor, 152mm	28	1
170-800-202	Balfour Retractor	27	1
170-800-241	Gelpi Point Retractor	27	1
170-800-231	Malleable Retractor Blade, 38mm	26	1
170-800-232	Malleable Retractor Blade, 51mm	26	1
170-800-233	Malleable Retractor Blade, 76mm	26	1
170-800-915	Case, 537x250x100mm	20	1

Blades, accessories and case


170-800-913 SET

Article No.	Description	Page	Pieces
170-800-100	Table Post for fixation at OP table	22	1
170-800-103	Horizontal Flex Bar	22	1
170-800-102	Horizontal Bar	22	1
170-800-032	Frame, oval, medium, segmented	24	1
170-800-012	Frame, oval, medium	24	1
170-800-914	Case, 537x250x197mm	20	1


Table fixation and ring case




More cases on demand.
Please contact yor area manager.



170-800-097	200mm, 8"
170-800-100	300mm, 12"
170-800-099	460mm, 18"
170-800-098	610mm, 24"



170-800-102	438mm, 17 1/4 "
-------------	-----------------




170-800-103	470mm, 18 1/2 "
-------------	-----------------




170-800-096	730mm, 28 3/4 "
-------------	-----------------



170-800-101	305mm, 12"
170-800-095	460mm, 18"




170-798-046	365mm, 14 2/5 "
-------------	-----------------




170-800-104	25x102mm, 1x4"
-------------	----------------



Fiber optic Clip	170-798-072
------------------	-------------




Fiber Optic Clip	170-800-071
------------------	-------------



Fiber optic Clip	170-800-070
------------------	-------------




Flexible Fiber Optic	170-798-060
----------------------	-------------




Ratchet, rotatable	170-800-106
--------------------	-------------



Standard Ratchet Mechanism	170-800-113
----------------------------	-------------



IMA Ratchet Mechanism	170-798-047
-----------------------	-------------



Ratchet Extender	170-800-105
------------------	-------------




Tilt Ratchet Mechanism*	170-800-110
-------------------------	-------------



Self-Retaining Tilt Ratchet*	170-800-111
------------------------------	-------------




Rotate- and tiltable Ratchet Mechanism*	170-800-112
---	-------------



Self-retaining rotate- and tiltable Ratchet Mechanism*	170-800-114
--	-------------

* detailed ratchet function illustrated on Page 28-29



Frame, round

170-800-021	small	∅ 248mm, 9 3/4 "
170-800-022	medium	∅ 292mm, 11 1/2 "
170-800-023	large	∅ 363mm, 14 1/4 "



Ring Segment, round, 2 pieces

170-800-041	small	∅ 254mm, 10"
170-800-042	medium	∅ 300mm, 12"




Frame, oval

170-800-011	small	340 x 245mm, 13 1/2 x 9 1/2 "
170-800-012	medium	410 x 294mm, 16 x 11 1/2 "
170-800-013	large	500 x 340mm, 19 3/5 x 13 1/2 "




Frame, oval, segmented

170-800-031	small	393 x 254mm, 15 1/2 x 10 "
170-800-032	medium	538 x 300mm, 21 x 12 "
170-800-033	large	544 x 354mm, 21 1/2 x 14 "




Elliptical Frame, segmented

170-800-035	369 x 381mm, 14 1/2 x 15 "
-------------	----------------------------




Vaginal Frame, oval, 20° angled

170-800-051	203 x 267mm, 8 x 10 1/2 "
-------------	---------------------------




Vaginal Frame, oval, segmented

170-800-050	305 x 241mm, 12 x 9 1/2 "
-------------	---------------------------



Ring Segment, 2 pieces

170-800-061	small	100mm, 4"
170-800-062	medium	150mm, 6"
170-800-063	large	200mm, 8"




Malleable Retractor Blade, with lip, Titan coated



with lip



Malleable Retractor Blade, with lip




Extra malleable Retractor Blade, with lip



with lip



Malleable Retractor Blade



Titanium Retractor Blade


170-800-229	38 x 279mm, 1 1/2 x 11"
170-800-230	51 x 279mm, 2 x 11"
170-800-227	76 x 279mm, 3 x 11"
170-800-228	102 x 279mm, 4 x 11"

170-800-237	25 x 102mm, 1 x 4"
170-800-238	38 x 168mm, 1 1/2 x 6 3/4"
170-800-239	51 x 168mm, 2 x 6 3/4"
170-800-240	76 x 168mm, 3 x 6 3/4"

170-800-247	25 x 102mm, 1 x 4"
170-800-248	38 x 168mm, 1 1/2 x 6 3/4"
170-800-249	51 x 168mm, 2 x 6 3/4"
170-800-250	76 x 168mm, 3 x 6 3/4"

170-800-231	38 x 152mm, 1 1/2 x 6"
170-800-232	51 x 152mm, 2 x 6"
170-800-233	76 x 152mm, 3 x 6"
170-800-234	76 x 254mm, 3 x 10"
170-800-235	102 x 215mm, 4 x 8 1/2"
170-800-236	51 x 254mm, 2 x 10"

170-800-360	25 x 102mm, 1 x 4"
170-800-361	25 x 127mm, 1 x 5"
170-800-362	25 x 152mm, 1 x 6"
170-800-363	25 x 178mm, 1 x 7"
170-800-364	25 x 203mm, 1 x 8"
170-800-365	51 x 102mm, 2 x 4"
170-800-366	51 x 127mm, 2 x 5"
170-800-367	51 x 152mm, 2 x 6"
170-800-368	51 x 178mm, 2 x 7"
170-800-369	51 x 203mm, 2 x 8"




Gelpi Point Retractor



Scoville Hook




Rake Retractor



IMA Retraktor



Rectal Retractor Blade



Balfour Retractor Blade



Hook Retraktor Blade

170-800-241	165mm
-------------	-------

170-800-321	small 25 x 40mm
170-800-322	large 25 x 80mm


170-800-343	3 prongs, 19 x 17mm
170-800-344	4 prongs, 9 x 19mm
170-800-345	6 prongs, 15 x 29mm

170-798-048	330mm, 13"
-------------	------------


170-800-311	small 22 x 72mm, 3/4 x 2 3/4"
170-800-312	large 23 x 98mm, 3/4 x 3 3/4"

170-800-202	75 x 75mm, 3 x 3"
170-800-201	89 x 70mm, 3 1/2 x 2 3/4"

170-800-251	26 x 25.5mm, 1 x 1"
170-800-252	51 x 25.5mm, 2 x 1"
170-800-253	102 x 25.5mm, 4 x 1"



170-800-221	64 x 143mm, 2 1/2 x 5 5/8"
-------------	----------------------------



170-800-222	90 x 251mm, 3 1/2 x 10"
-------------	-------------------------



170-800-208	25 x 38mm, 1 x 1 1/2"
170-800-209	25 x 63mm, 1 x 2 1/2"
170-800-210	51 x 38mm, 2 x 1 1/2"
170-800-211	51 x 63mm, 2 x 2 1/2"
170-800-212	51 x 51mm, 2 x 2"
170-800-214	51 x 76mm, 2 x 3"
170-800-215	51 x 102mm, 2 x 4"
170-800-216	51 x 127mm, 2 x 5"
170-800-217	51 x 152mm, 2 x 6"
170-800-218	51 x 200mm, 2 x 8"



170-800-351	25 x 76mm, 1 x 3"
170-800-352	25 x 102mm, 1 x 4"
170-800-353	25 x 127mm, 1 x 5"
170-800-354	25 x 152mm, 1 x 6"
170-800-355	25 x 180mm, 1 x 7"
170-800-356	25 x 203mm, 1 x 8"




170-800-291	19 x 25mm, 3/4 x 1"
170-800-292	19 x 51mm, 3/4 x 2"
170-800-293	25 x 32mm, 1 x 1 1/4"



170-800-261	small, 25 x 89mm, 1 x 3 1/2"
170-800-262	large, 25 x 127mm, 1 x 5"




170-800-331	small, 60 x 17mm, 2 1/3 x 2/3"
170-800-332	medium, 85 x 25mm, 3 1/3 x 1"
170-800-333	large, 95 x 50mm, 3 3/4 x 2"



170-800-281	38 x 51mm, 1 1/2 x 2"
170-800-282	40 x 25mm, 1 1/2 x 1"
170-800-283	65 x 25mm, 2 1/2 x 1"
170-800-284	122 x 20mm, 2 3/4 x 3/4"
170-800-285	122 x 25mm, 2 3/4 x 1"



170-800-271	25 x 76mm, 1 x 3"
170-800-272	25 x 114mm, 1 x 4 1/2"



1



- self retaining
- Retractor rotatable
- Retractor tiltable




170-800-114

2



- self retaining
- Retractor tiltable




170-800-111

3



- Retractor rotatable
- Retractor tiltable




170-800-112

4



- Retractor tiltable



170-800-110

Media Code Number: 021241

© 10/2019 Copyright by RZ Medizintechnik GmbH · Unter Hasslen 20 · D-78532 Tuttlingen. All rights reserved.
Modifications of technical specification or design are subject to change without prior notice.