

Welch Allyn Pocket Junior Otoscopes and Ophthalmoscopes

Small Scopes That Make a Big Difference

Welch Allyn Pocket Junior diagnostic instruments are valuable diagnostic tools for any office practice. Built with the same technology as Welch Allyn's other high-quality products, provide the illumination needed to help you diagnose at an affordable price. And with its attractive styling and convenient size, it will quickly become a big part of your practice.



Welch Allyn Pocket Junior Otoscope

- Lightweight and compact
- Halogen illumination provides brighter, whiter, longer lasting light
- Fiber optic light transmission
- Wide angle magnifying lens
- Insufflation port for pneumatic otoscopy
- Standard AA batteries (not included) work with interchangeable heads for otoscope and ophthalmoscope



Welch Allyn Pocket Junior Ophthalmoscope

- Lightweight and compact
- Halogen illumination provides brighter, whiter, longer lasting light
- Sealed optics keep instrument clean and dust free
- Aperture dial offers six selections
- Glass condensing lens
- Includes red-free and cobalt blue filter
- Standard AA batteries (not included) with interchangeable heads for otoscope and ophthalmoscope



Welch Allyn Pockot Junior Sot

Pocket Junior Set — A Winning Combination

Put the power, convenience, and versatility of a Welch Allyn Pocket Junior Set including an otoscope and ophthalmoscope to work for you. This set will provide the valuable diagnostic tools you need at an affordable price, and it comes in an attractive, yet durable case that's easy to carry.

Ordering Information



Welch Allyn Pocket Junior Otoscope

- 22840 Pocket Junior Otoscope complete with AA handle*
- 22841 Pocket Junior Otoscope complete with AA handle* and case
- 03400 Replacement Lamp
- 21504 Otoscope Insufflator Bulb
- 52432 Disposable Specula 2.5 mm
- 52434 Disposable Specula 4 mm

Welch Allyn Pocket Junior Ophthalmoscope

12850 Pocket Junior Ophthalmoscope complete with AA Handle*
12851 Pocket Junior Ophthalmoscope with AA Handle* and Case





Welch Allyn Pocket Junior Set

95001 Pocket Junior Complete Set: Otoscope with 25 specula, Ophthalmoscope, 2 Handles and Case

Welch/Allyn

*Batteries not included

Specifications: Welch Allyn Pocket Junior Otoscopes and Ophthalmoscopes

	Pocket Junior Otoscope	Pocket Junior Ophthalmoscope
Illumination	Halogen	Halogen
Transmission	Fiber Optic to distal tip	_
Batteries	2 AA (not included)	2 AA (not included)
Viewing Lens	Wide angle and removable	Glass Condensing Lens
Weight with Batteries	132 g (4.65 oz)	129.5 g (4.56 oz)
Height	16.33 cm (6.43 in)	16.02 cm (6.31 in)
Diopter Range	—	-25 to +22
Handle Material	All metal construction	All metal construction
Insufflation Port	Yes	—
Aperture Positions	—	6 including Slit and Fixation
Sealed	—	Yes
Filters	—	Cobalt Blue and Red-Free

