

GET250L, GET350L

Fluid Warming Cabinets

Single chamber Fluid warming cabinets. Interior surface temperatures are monitored and balance heat throughout the cabinets. The fluid chamber can be programmed to warm either irrigation fluids (IRR) or injection fluids (INJ) with separate temperature ranges provided depending on the choice selected. The fluid portion can hold bottles, bags or a combination of both.

IRR: 32 °C to 66 °C/ 90 °F to 150 °F)

INJ: 32 °C to 40 °C/ 90 °F to 104 °F)



Application

The Getinge GET250L and GET350L are designed to warm injection or irrigation fluids. This product is not to be used to warm blood or blood products.

Construction

- Stainless Steel exterior and white epoxy-coated interior
- Fully gasketed windowed doors with magnetic closure
- Four non-skid rubber feet are standard
- Right-hinged standard
- Shelf Standard on GET350L, no shelf on GET250L

Dimensions

Model	Exterior (H x W x D)	Interior (H x W x D)	Net Weight
GET250L	22"(557mm) x 18.5"(470mm) x 23.5"(597mm)	15.25"(387mm) x 15.5"(394mm) x 18.75"(476mm)	67 lbs (30 kg)
GET350L	28"(709mm) x 18.5(470mm) x 23.5"(597mm)	21.25"(540mm) x 15.5"(394mm) x 18.75"(476mm)	91 lbs (41 kg)

Capacity

Model	Bottles	Bags	Blankets
GET250L	12	16	N/A
GET350L	24	24	N/A

Installation

Stationary or mobile warming cabinets

Number of doors

- Single door with window

Installation alternatives

- Bumper and Caster Assembly
- Caster and Plate Assembly
- Solid Door
- Left-hinged
- One or two baskets
- Mobile Cart

Other installation options

- Manual lock
- Seismic Mounting brackets (specific installation requirements apply)
- Stacking Kit

Clearance Requirements

- 3" (76mm) from rear
- 1" (25mm) from top and sides
- 3/4" (19mm) from bottom

Safety features

- The controller stores all settings and resumes operation using these settings in the event of a power failure.
- A warming shut-off system, separate from the controller, prevents overheating.

Other features

- Access point at rear panel allows the addition of data logging or temperature management hardware.
- WarmSafe™ heating system incorporates a multiple zone warming technology (Patent No: US 8,217,316; US 8,581,152) that heats where and when it is needed.
- Operate in Celsius or Fahrenheit
- LED display
- Built-in speaker for audible feedback
- LED interior lighting

Getinge is a registered trademark of Getinge AB, its subsidiaries, or affiliates in the United States or other countries • Copyright 2019
Getinge AB or its subsidiaries or affiliates • All rights reserved • MCV00076066 REV B • EXP 02/21

Sales Office: Getinge USA Sales LLC, 1777 East Henrietta Road, Rochester, NY 14623

Manufactured By: Enthemics Medical Systems, W164 N9221 Water Street, Menomonee Falls, WI 53051